



CALL-DEX

MOBILE PHONES COMPATIBILITY GUIDE



By choosing Widex hearing aids, you are choosing products from a company that has been WindMade certified. WindMade is the first global consumer label identifying companies that use wind power.

WWW.WIDEX.COM







Widex and DEX are trademarks of Widex A/S.



2015-05
9 502 3815 001 #01

Widex A/S, Nymoellevvej 6, DK-3540 Lyngø



NOTE: If you use the phone for conversations for more than 3 hours a day, we recommend changing the battery once a week. If you use it less than 3 hours a day, we recommend changing the battery once every second week.

CALL-DEX uses a size 10 battery.
Battery is not included in the case.



WHICH MOBILE PHONE CAN I USE WITH CALL-DEX?

CALL-DEX is compatible with mobile phones from a wide range of manufacturers including Apple, HTC, Sony, Samsung and Nokia. The list on the right names the specific phone models.



ACER

- Cloud Mobile
- Liquid E3
- Liquid Gallant
- Liquid Glow

ALCATEL

- OT Scribe HD
- OT-916
- OT-918D
- OT-T*POP

APPLE

- iPhone 4
- iPhone 4s
- iPhone 5
- iPhone 5c
- iPhone 5s
- iPhone 6
- iPhone 6 Plus

BLACKBERRY

- Classic
- Q10
- Q5
- Z10

BQ

- Aquaris E4.5

DORO

- Liberto 810
- Liberto 820

EMPORIA

- Life

HTC

- Desire 601
- Desire 816
- Desire C
- Desire HD
- Incredible S
- One
- One (M8)
- One M9
- One Mini 2
- One S
- One SV
- One V
- One X
- Sensation XE
- Windows Phone 8S
- Windows Phone 8X

HUAWEI

- Ascend G510
- Ascend Mate
- P7
- Ascend Y300
- Ascend Y530
- Honor 6

LG

- G Flex
- G2
- G2 Mini
- G3
- G3S
- Nexus 4
- Nexus 5
- Optimus 2X
- Optimus 4X HD
- Optimus L3
- Optimus L7 II
- Optimus Me
- Optimus One

MICROSOFT

- Lumia 535

MOTOROLA

- MILESTONE 2
- Moto E
- Moto G
- Moto X
- MOTOLUXE
- Nexus 6
- RAZR
- RAZR i

NOKIA

- 113
- 215
- 700
- Asha 202
- Asha 303
- Asha 503
- Lumia 1020
- Lumia 1320
- Lumia 1520
- Lumia 520
- Lumia 530
- Lumia 610
- Lumia 620
- Lumia 625
- Lumia 630
- Lumia 735
- Lumia 820
- Lumia 920
- Lumia 925
- Lumia 930
- X

SAMSUNG

- ATIV S
- Galaxy A3
- Galaxy Ace
- Galaxy Ace 2
- Galaxy Ace Plus
- Galaxy Ace Style
- Galaxy Alpha
- Galaxy Fame
- Galaxy Fit
- Galaxy Gio
- Galaxy Mega 6.3
- Galaxy Nexus
- Galaxy Note
- Galaxy Note 3
- Galaxy Note 4
- Galaxy Note Edge
- Galaxy Note II
- Galaxy Pocket Neo
- Galaxy S6
- Galaxy S II
- Galaxy S III
- Galaxy S III Mini
- Galaxy S Plus
- Galaxy S4
- Galaxy S4 Mini
- Galaxy S5
- Galaxy S 5 Mini
- Galaxy W
- Galaxy Xcover

SEMC

- Xperia Arc
- Xperia Arc S
- Xperia Mini Pro
- Xperia X8

SONY

- Xperia E
- Xperia E1
- Xperia E3
- Xperia go
- Xperia L
- Xperia P
- Xperia SP
- Xperia T
- Xperia T3
- Xperia Tipo
- Xperia U
- Xperia Z
- Xperia Z1
- Xperia Z1 Compact
- Xperia Z2
- Xperia Z3 Compact

WILKO

- Cink Peax
- Wax

ZTE

- Blade III
- Blade Q Mini
- Grand X
- Kis
- Libra

Disclaimer: All mobile phones in this list have been tested comprehensively by Widex, and have been found to fulfil the minimum requirements for CALL-DEX compatibility at the time of testing. However, Widex takes no responsibility for such compatibility in general. Specifically, Widex shall not be responsible for incompatibility resulting from future updates of the mobile phones, nor from events of incompatibility not covered by and/or found in these tests. Widex takes no financial responsibility in terms of reimbursements of mobile phones which do not work as expected.